



Printed or electronic versions of this document not accessed directly from the designated Hillrom Controlled Quality Information System are For Reference Only.

## ***Customer Service Bulletin***

<b>Product:</b>	SmartCare™ Remote Management (SCRM)		
<b>Subject:</b>	Introducing Hillrom's <i>SmartCare™ Remote Management</i> (SCRM replaces the Welch Allyn Service Hub)		
<b>HW Version(s) Affected:</b>	N/A	<b>SW Version(s) Affected:</b>	6.2.1

<b>Classification:</b>	Information Only		
<b>Distribution:</b>	<input checked="" type="checkbox"/> Customer Care	<input checked="" type="checkbox"/> Product Service	<input checked="" type="checkbox"/> Field Service
	<input type="checkbox"/> ASPs <input type="checkbox"/> Distributors	<input checked="" type="checkbox"/> Customers	<input type="checkbox"/> Company Confidential

### **Summary:**

#### **Introducing Hillrom's *SmartCare™ Remote Management***

SmartCare Remote Management helps our customers service assets remotely. It supports the following service-related tasks:

- Remote update configuration
- Remote upgrade asset firmware
- Remote asset log retrieval
- Service data Reports
- Review asset preventive maintenance (PM) due list
- Remote asset location tracking

Existing Welch Allyn Service Hub users will be automatically redirected to the new SmartCare Remote management URL:

<https://smartcareremotemanagement.hillrom.com>



## Version 6.2.1 Details:

- Added Progressa bed support to SCRM
- Single Sign-on (SSO) option for customer login to Smart Care Remote Management  
The single sign-on method allows customers to authenticate their credentials through the company's existing Active Directory (Windows login). This method eliminates the need to provide a username and password at every login. Contact your Hillrom representative for setup details.  
**Single sign-on (customer option) requirements:**  
For use with cloud-based Microsoft Azure Active Directory only (no on-prem legacy). The user authentication must be compatible with the SmartCare Remote Management single sign-on option.
- Update Activity report to remove all the fields and add the columns as specified in the details (Report columns will differ based upon asset name)
- Bug fix: When editing APLocation details, Campus and Building information swapped.
- Bug fix: Customer Portal: duplicate records in csv download.

**System Compatibility:** Supported Operating Systems for PartnerConnect Agent (applicable assets only. See Admin Guide for more detail):

- Windows 10 (32 and 64-bit)
- Windows Server 2016 (64-bit)
- Windows Server 2019 (64-bit)

## Web Browser Compatibility:

*Note: The preferred browser is Google Chrome.*

- Google Chrome 92 or higher
- Microsoft Edge 92 or higher
- Apple Safari iOS 14

**Supported Hillrom Applications:**

	Version
PartnerConnect Agent	3.0.0 or later
Welch Allyn Service Monitor	All
Connex CS Server	All except 1.8.2

**Supported Hillrom Devices:**

	SCRM
Connex Spot Monitor	1.24 or later
CVSM/CIWS	2.x or later except version 2.40.x
RV700	1.30.00-A0002
Centrella Bed	1.36.x or later
Progressa Bed	22021501

See SmartCare Remote management URL <https://smartcareremotemanagement.hillrom.com>

for the latest Administration Guide, User Guide, and PartnerConnect Installer.

**END OF BULLETIN**

Version	Sec, Pg, Para Changed	Change Made	Date Version Created	Version Created By (initials)
A	N/A	Initial Release	2021-06-17	DCS
B	All	Changes applied for SCRM. 6.1.1 updates	2021-10-13	DCS
C	All	Changes applied for SCRM. 6.2.1 updates	2022-03-01	DCS